
# MRI

---

Manuale Tecnico MRI  
MRI Technical Manual  
Manuel technique IRM  
Manual técnico MRI

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## CONTENTS

<b>1.</b>	<b>INTRODUCTION</b> .....	<b>9</b>
<b>2.</b>	<b>MEDICO S.P.A. MR-CONDITIONAL SYSTEM CONDITIONS FOR USE</b> ....	<b>9</b>
<b>3.</b>	<b>GENERAL INFORMATION ON MRI MACHINES</b> .....	<b>10</b>
3.1	TYPES OF ELECTROMAGNETIC FIELDS GENERATED BY MRI MACHINES ....	10
3.2	POTENTIAL INTERACTIONS FOR IMPLANTED PACING SYSTEMS IN THE MRI ENVIRONMENT .....	10
<b>4.</b>	<b>MEDICO MRI-CONDITIONAL PACING SYSTEM CONTRAINDICATIONS</b> .....	<b>10</b>
<b>5.</b>	<b>POTENTIAL ADVERSE EVENTS CORRELATED TO THE USE OF MEDICO MR-CONDITIONAL PACING SYSTEM</b> .....	<b>11</b>
<b>6.</b>	<b>CARDIOLOGY-SPECIFIC WARNINGS AND PRECAUTIONS RELATED TO THE USE OF MEDICO MR-CONDITIONAL PACING SYSTEM</b> .....	<b>11</b>
6.1	CARDIAC EQUIPMENT REQUIREMENTS .....	11
6.2	PATIENT INFORMATION .....	11
<b>7.</b>	<b>RADIOLOGY-SPECIFIC WARNINGS AND PRECAUTIONS RELATED TO MEDICO MR-CONDITIONAL PACING SYSTEM</b> .....	<b>12</b>
7.1	REQUIREMENTS OF THE MR EQUIPMENT.....	12
7.2	MRI SCAN PROTOCOL REQUIREMENTS .....	12
7.3	PATIENT SAFETY DURING THE MRI SCAN .....	12
7.4	IMAGE QUALITY .....	12
<b>8.</b>	<b>OPERATIONS TO BE DONE BEFORE MAGNETIC RESONANCE SCAN</b> ... <b>12</b>	
8.1	REQUIRED PATIENT CARE .....	12
8.2	IDENTIFICATION OF MEDICO MRI PACING SYSTEM COMPONENTS.....	12
8.3	PREPARING A PATIENT WITH A MEDICO MRI PACING SYSTEM FOR AN MRI SCAN .....	12
<b>9.</b>	<b>INITIATING "MRI MODALITY" FEATURE</b> .....	<b>13</b>
9.1	MEDICO MRI PACING SYSTEM INTEGRITY VERIFICATION.....	13
9.2	"MRI MODALITY" FEATURE ACTIVATION .....	14
<b>10.</b>	<b>DURING "MRI MODALITY" OPERATION</b> .....	<b>14</b>
10.1	DEVICE CONSIDERATIONS.....	14
10.2	RADIOLOGY CONSIDERATIONS DURING THE MRI SCAN .....	14
<b>11.</b>	<b>AFTER THE MRI SCAN</b> .....	<b>15</b>
11.1	RESTORING THE DEVICE TO THE PRE-MRI CONFIGURATION.....	15

## 1. INTRODUCTION

The Eos DR MRI and Eos D MRI pulse generators manufactured by Medico S.p.A. are designed to allow patients to be safely scanned by an MRI machine (MRI exam), when used according to the specific MRI conditions for use, reported in this document. Such conditional security is defined by standardized terminology (standard ASTM F2503 2013) as *MR-conditional*.

Before conducting an MRI scan on a patient implanted with a Medico MRI certified pacing system, it is necessary to read carefully in this manual information about procedure, warnings and specific Magnetic Resonance precautions.

Possible doubts or questions about the subjects contained in this manual can be clarified contacting Medico S.p.A.

Refer to the appropriate Medico S.p.A. pacemaker or lead user's manual for non-MRI related instructions for use.

**"MRI-CONDITIONAL" qualification does not apply to Eos DR MRI or Eos D MRI pacemakers in themselves, but to a system necessarily comprising Eos DR MRI or Eos D MRI pacemaker connected to a 400 MRI (ventricular passive fixation bipolar pacing lead) and a 400J MRI (atrial "J" shaped passive fixation bipolar pacing lead).**

## 2. MEDICO S.p.A. MR-CONDITIONAL SYSTEM CONDITIONS for USE

Medico S.p.A. has performed accurate and extensive tests, which have demonstrated that its own pacing system can be used safely in the MRI environment provided that the instructions reported in this manual are appropriately followed.

### Cardiology checklist:

- No other active or abandoned implanted devices, components, and/or accessories present in the patient.
- The implanted system must be composed solely by a Medico MRI certified system (EOS DR MRI or EOS D MRI pacemaker with 400 MRI and 400J MRI leads).  
**WARNING:** Do not bring any system components that are not marked MR-safe  or MR-conditional  in the magnet room.
- Pacing system must be set to "MRI Modality"<sup>2</sup> during the scan.
- Pacing system implanted for more than 6 weeks.
- Pulse generator implant location restricted to left or right sub-clavicular pectoral region.
- Pacing threshold in bipolar  $\leq 2.0$  V at a pulse width of 0.43 ms.
- No damaged lead or compromised system integrity. (Note: the lead impedance value must be  $\geq 200$  ohms ( $\Omega$ ) and  $\leq 2000$  ohms ( $\Omega$ )).

<sup>2</sup> During MRI Modality, pacing mode is temporarily set to VOO or OOO (bipolar), the ventricular output energy is increased and pacing rate is set among 80 to 100 min<sup>-1</sup>. The MRI Modality will last for a programmable time. When the MRI Modality is terminated, either by the programmer or because this time has expired, the pacemaker will revert to previous conditions **except for the ventricular output energy**, that will remain high until it is manually reprogrammed.

- Pacemaker battery impedance lower than 5000 ohms ( $\Omega$ ).

### **Radiology checklist:**

Patients implanted with a Medico MR-Conditional pacing system can be scanned under the following conditions:

- **The MRI scanning is allowed for head, extremities and lumbosacral region. The thorax is therefore excluded.**
- Transmit-receive or transmit-only head and local RF coils that do not require placement on the thorax.  
**WARNING:** Do not use transmit-receive or transmit-only local RF coils positioned on the thorax.  
**WARNING:** Do not use transmit-receive or transmit-only RF body coils.
- Horizontal cylindrical closed-bore magnet is used.
- An MRI scanner with a static magnetic field of 1.5 T (tesla) must be used.
- Spatial gradient of the static magnetic field must not exceed 20 T/m or 2000 G/cm.
- Gradient systems must not exceed a slew rate of 200 T/m/s per axis.
- The MRI scanner Operating Mode must be: "Normal Operating Mode".
- Patient must be scanned only in supine position.
- The distance between the implant and the wall of the bore must be at least 5 cm.
- Patient must not have elevated body temperature or compromised thermoregulation at time of scan.  
**WARNING:** Do not scan a patient with an elevated body temperature.
- Continuous patient monitoring must be provided. Maintain visual and verbal contact with the patient, and continuous hearth rate monitoring by pulse oximetry or electrocardiography (ECG).
- In non-clinical testing, the maximum image artifact was seen on the gradient echo pulse sequence and extended from the boundary of the implant components by:

72 mm	from the	Pulse Generator	EOS D, DR MRI
11 mm	from the	Atrial Lead	400J MRI
13 mm	from the	Ventricular Lead	400 MRI

- The total duration of the Magnetic Resonance exam must not exceed 30 minutes.

## **3. GENERAL INFORMATION ON MRI MACHINES**

### **3.1 Types of electromagnetic fields generated by MRI machines**

An MRI machine typically produces three types of electromagnetic fields which may interact with implanted device systems. All of these fields are essential to produce an MRI image. These fields are defined as follows:

Static magnetic field - Steady state non-varying magnetic field that is usually always present around an MRI machine, even when no diagnostic exam is underway

Gradient magnetic fields - Low-frequency pulsed magnetic fields that are present only during a diagnostic exam. Three mutually orthogonal gradient magnetic fields are used by MRI equipment to form the image.

RF field - Pulsed radio frequency (RF) field that is present only during a diagnostic exam. This kind of field can be produced by different transmission RF coils, such as a whole body transmit coil or a local coil (for example, a transmit head coil or an extremity coil). Generation of the RF field with the whole-body coil is not allowed with the Medico S.p.A. MR-conditional pacing system.

### **3.2 Potential interactions for implanted pacing systems in the MRI environment**

Medico MR-conditional pacing system is designed to mitigate the potential interactions described in this section.

Magnetic field interactions - The electrically conductive material of an implanted system may apply force, vibration, and torque effects (moment) as a result of the static magnetic field and gradient magnetic fields produced by an MRI scanner. The pacing system under discussion is designed to mitigate these effects and reduce the mechanical stress on the implanted system and tissue interface. Patients may feel a vibration or a slight movement sensation near the site of the device implant, inside or in proximity of the MRI scanner.

Induced stimulation - The RF and gradient magnetic fields produced by an MRI scanner generate energy that could cause potentially unintended heart stimulation. The Medico pacing system is built to limit impulse voltage and pulse width induced on the leads, to mitigate the risk of cardiac capture.

Electrodes and case heating - The RF fields generated by an MRI scanner induce voltages that may cause heating at the lead electrodes. This warming could damage tissue that surrounds the electrodes, thus compromising pacing and sensing thresholds at that site. By excluding the use of the RF field generated by the whole body coil and not allowing local coils that are not placed around the thorax, Medico MR-conditional leads are not affected by the heating due to the RF field. An increase of the temperature of the pacemaker case could occur, but test have demonstrated that dangerous temperatures can't be reached.

Image artifacts and distortion - Medico MR-conditional leads have demonstrated to grant minimal image distortion for contiguous areas surrounding the implanted leads, when the pacemaker is out of the field of view. The presence of the device near the field of view could determine a significant image distortion. When selecting the field of view and image-acquisition parameters, image artifacts and distortion (resulting from the presence of the device and the leads near the field of view) must be considered.

#### 4. MEDICO MRI-CONDITIONAL PACING SYSTEM CONTRAINDICATIONS

- MRI scan is contraindicated for patients with previously not MR-conditional implanted medical devices even if inactive, partially removed or abandoned.
- MRI scan is contraindicated for patients with damaged, not connected, abandoned, intermittent or partially removed leads. Leads originally part of an MR-conditional system but disconnected/abandoned at a later stage may cause strong harm to the patient, so they must be considered as a contraindication source. MRI scan is strongly contraindicated for patients with failed/compromised leads or leads or with intermittent functionality.
- MRI scan is contraindicated for patients with a Medico MRI pacing system implanted for less than 6 weeks.
- MRI scan is contraindicated for patients with a Medico MRI pacing system implanted in locations not restricted to the left and right sub-clavicular pectoral region.
- MRI scan is contraindicated for patients who do not have a complete Medico MRI pacing system which must include a Eos DR MRI or Eos D MRI pacemaker with atrial 400J MRI and ventricular 400 MRI leads.
- MRI scan is contraindicated for patients with bipolar ventricular pacing capture threshold values  $\geq 2.0$  V at a pulse width of 0.43 ms.  
**Note:** If all other MRI scan pre-requirements are satisfied, patients experiencing atrial fibrillation may undergo an MRI examination.
- MRI scan is contraindicated whenever it is not possible to have a patient monitoring system (ECG, pulse oximetry, non-invasive blood pressure measurement) available or whenever it is not possible to maintain visual and verbal contact with the patient.
- MRI scan is contraindicated for patients that are affected by diaphragmatic stimulation when the device is programmed to pace in MR-conditional mode, i.e. to asynchronous pacing with an output amplitude of 5.0 V and a pulse width of 0.98 ms. OOO pacing mode can be used in this case, if appropriate to the patient.
- MRI scan is contraindicated for patients with a lead impedance value  $\leq 200$   $\Omega$  or  $\geq 2000$   $\Omega$ .
- MRI scan is contraindicated with patient positioned differently from the supine position; inside the MRI bore, a patient with an implanted Medico MR-conditional pacing system must not be positioned prone or on his or her side (lateral decubitus position).
- The use of the whole body coil for transmission is contraindicated, transmission can only be from a local coil that does not include the thorax.
- MRI scan is contraindicated for patients who do not tolerate the OOO mode (stimulation switch off) or the asynchronous ventricular stimulation. These are the only pacing modes that are applied by the EOS DR MRI and EOS D MRI pacemakers during the MR-conditional mode.
- MRI scan is contraindicated in patients affected by body temperature higher than normal (fever) or with poor control of the body temperature.

#### 5. POTENTIAL ADVERSE EVENTS CORRELATED TO THE USE OF MEDICO MR-CONDITIONAL PACING SYSTEM

Medico MR-conditional pacing system and its conditions for use are designed to minimize any potential adverse event that may cause harm to the patient. In the MRI environment the following potential adverse events may occur:

- Lead electrodes heating and tissue damage, causing loss of sensing and/or capture.
- Device heating causing tissue damage in the implant pocket and/or patient discomfort/ache.
- Induced currents on leads causing continuous capture and/or VT/VF or hemodynamic collapse; damaged device or leads, resulting the system to fail to detect or treat irregular heartbeats or causing the system to treat patient's condition incorrectly.
- Damage of the functionality or the mechanical integrity of the device, causing the inability of the device to communicate with the programmer.
- Vibration or movement of the device or leads, resulting in device or lead dislodgment.
- Competitive pacing and potential VT/VF induction due to asynchronous pacing when the MRI modality is enabled.

#### 6. CARDIOLOGY-SPECIFIC WARNINGS AND PRECAUTIONS RELATED TO THE USE OF MEDICO MR-CONDITIONAL PACING SYSTEM

##### 6.1 Cardiac equipment requirements

"MRI Modality" configuration - Security safeguard of the patient during each Magnetic Resonance session requires that the implanted pacing system exclusively consists of an Eos DR MRI or Eos D MRI pacemaker with 400 MRI and 400J MRI, respectively ventricular and atrial leads. A properly designed configuration has to be programmed preventively, that is accessible from the pacemaker management software installed in Olympia programmer. Please refer to the programmer user's manual for the operative details.

**Important note:** be aware that some patients may be susceptible to cardiac arrhythmia induced by competitive pacing, when an asynchronous pacing mode is selected. It is really important for such patients to first select a pacing rate that avoids competitive pacing and then shorten as possible the duration of the asynchronous pacing.

##### 6.2 Patient information

System information and records - All relevant information about the components of the implanted Medico MRI pacing system, such as model names, serial numbers, and leads lengths must be recorded in the patient's record and set up on the programmer patient information screen. If available, the MRI patient record sticker has to be placed in the patient record too.

**Note:** it is important to confirm that the lead model and its length are recorded in the patient records.

**Patient ID card requirements** – All the patients with an implanted Medico MRI pacing system must carry a proper identification document, such as patient's ID card. This document must indicate that the patient is implanted with a Medico MRI pacing system.

## 7. RADIOLOGY-SPECIFIC WARNINGS AND PRECAUTIONS RELATED TO MEDICO MR-CONDITIONAL PACING SYSTEM

### 7.1 Requirements of the MR equipment

Verify that the MRI equipment requirements listed in this section are carefully satisfied during all MRI scans performed on patients with a Medico MR-conditional pacing system. In case of doubt or uncertainty regarding the features of the MRI machine used, please contact the MRI manufacturer or its representative.

**Acceptable MRI equipment operating features** – Safety and reliability of Medico MR-conditional pacing system have been tested for diagnostic exams using MRI equipment provided with the following operating characteristics:

- Hydrogen proton magnetic resonance imaging equipment with a static magnetic field of 1.5 T.
- RF excitation frequency in the amount of approximately 64 MHz.
- MRI systems for clinical use, with horizontal cylindrical closed-bore magnet.
- A gradient slew rate performance  $\leq 200$  T/m/s or less per axis.

### 7.2 MRI scan protocol requirements

The Magnetic Resonance must be used under "normal operating mode"; the "first level controlled" and "second level controlled" modes must not be used.

### 7.3 Patient safety during the MRI scan

**Defibrillator** – During the MRI scan an external defibrillator must be available.

**Patient monitoring** – During the MRI scan, a proper patient monitoring must be guaranteed. The Magnetic Resonance environment may interfere with the patient monitoring system, therefore the use of more than one of the following systems is strongly recommended:

- Electrocardiography;
- Pulse oximetry;
- Non-invasive blood pressure measurements.

In addition, visual and verbal contact with the patient must be maintained.

**Note:** If the patient's hemodynamic is compromised during the MRI scan, suspend the MRI examination and implement the proper measures to restore the patient's hemodynamic performance.

## 7.4 Image quality

Whenever the Medico MR-conditional pacing system is near the imaging area of interest, the quality of the images obtained with Magnetic Resonance could be deteriorated.

## 8. OPERATIONS TO BE DONE BEFORE MAGNETIC RESONANCE SCAN

### 8.1 Required patient care

Before performing a MRI scan, it is necessary to verify the following steps:

- Verify that the pacing capture threshold values (in bipolar) are  $\leq 2.0$  V at a pulse width of 0.43 ms.  
**Note:** Patients affected by atrial fibrillation may be scanned, provided all other pre-MRI scan requirements are satisfied.
- Verify the absence of diaphragmatic stimulation for patients whose device will be programmed to an asynchronous pacing mode when "MRI Modality" is on. The presence of diaphragmatic stimulation is a contraindication to MRI scan.

### 8.2 Identification of Medico MRI pacing system components

Eos DR MRI and Eos D MRI pacemakers are equipped with radiopaque MRI symbols. It is therefore recommended to verify whether the device is an Eos DR MRI or Eos D MRI model by taking an x-ray of the implanted system. An x-ray examination also allows to ascertain whether the patient has implantable devices that are contraindicated for an MRI scan.

The 400 MRI and 400J MRI leads are not provided with the radiopaque symbol. It is therefore essential to verify whether the leads are of these models in the patient record at the hospital, or in the documents in the hands of the patient. Refer possibly also to the pacemaker and leads user's manuals for further information.

**Note:** The patient records or patient ID cards (if present) must be accurate and complete. The patient could have other MRI contraindicated implantable devices that are not reported on the Medico patient ID card.

### 8.3 Preparing a patient with a Medico MRI pacing system for an MRI scan

Before performing an MRI scan on a patient with a Medico MRI pacing system, it is necessary to complete the following tasks:

- Verify that the pacing system has been implanted for more than 6 weeks: the 6-weeks post-implant waiting period is essential to obtain the device and lead fixation in the patient. Through a proper maturation of the lead-tissue interface, the stability of the capture threshold is indeed increased and the risk of unintended pacing due to interference generated by the MRI system is mitigated.
- Verify that the pacemaker was implanted in the pectoral region: implant locations different from left or right pectoral region have not been analysed, so safety can't be guaranteed.

- Verify that no other implanted devices are present: interactions with all other implanted devices (active or unused) have not been tested by Medico. MRI scan can be considered safe for the patient only if there are no other implanted devices.
- Verify that the leads are Medico MRI leads: different leads have not been tested by Medico, so they can not be considered safe. MRI scan can be considered safe for the patient if and only if the implanted system comprises an Eos DR MRI or Eos D MRI pacemaker connected to 400 MRI (ventricular passive fixation bipolar pacing lead) and 400J MRI (atrial "J" shaped passive fixation bipolar pacing lead).
- Verify that the leads are electrically intact: any kind of fracture or other damage to the leads could alter the electrical properties of the pacing system and the system will be therefore unsafe for an MRI scan. Patients with damaged leads may be harmed if an MRI scan is performed.
- Verify that abandoned or additional leads are not present: all the other sort of leads (different from the ones above) same as disused leads have not been tested by Medico. MRI scan can be considered safe for the patient only if no other implanted leads are present.
- Verify that no lead extenders or adaptors are present: lead extenders or adaptors are not considered part of the Medico MRI pacing system, and have not been tested for an MRI scan. Their influence on the electrical characteristics of the implant is unknown and could be significant.

**Maximum pacing capture threshold** – Before performing an MRI scan, the atrial and ventricular pacing capture thresholds of the patient must be evaluated. An MRI scan is **contraindicated** if either the atrial or ventricular capture threshold in bipolar is  $\geq 2.0$  V at a pulse width of 0.43 ms.

**Pacing therapy during MRI scans** – Each patient must be examined to determine whether or not a pacing support is required during the MRI scan. If pacing support is needed, the pacing mode of the "MRI Modality" function must be set to VOO. Otherwise, if pacing support is not required, the pacing mode of the "MRI Modality" function should be set to OOO.

**Note:** The physician has the task and responsibility to select the most appropriate mode of operation for the patient.

**Pacing rate** – The most appropriate pacing rate for patients who require a pacing support must be established before performing an MRI scan. Pacing rate has to be selected appropriately by the physician (values among 80 to 100  $\text{min}^{-1}$ ) to avoid competitive pacing when the "MRI Modality" feature is on.

## 9. INITIATING "MRI MODALITY" FEATURE

**WARNING:** Do not bring the Medico programmer and its accessories into the magnet room, as it is neither MR-conditional nor MR-safe.

To perform a safe MRI scan on a patient with an implanted Medico MRI pacing system, the "MRI Modality" feature must be enabled using the Olympia programmer provided with a software version 09.0.0 or above.

To program the MRI scan feature, the user has to select the appropriate parameters for the patient before starting the MRI scan.

**Note:** Pacing mode and rate (if pertinent) must be programmed at discretion of the physician. It is possible to program an asynchronous pacing mode (VOO) or OOO mode (i.e., no pacing), depending on whether the patient needs pacing support or not. When the "MRI Modality" is enabled, sensed events are ignored by the device, independently from the programmed mode. Programmable parameters of MRI Modality are listed in the table below.

Parameter	Programmable values
Maximum time under MRI mode (hours)	12/24/36/48
MRI mode	VOO/OOO
Basic rate ( $\text{min}^{-1}$ )	80/90/100
Ventricular pulse amplitude (V)	5/6
Ventricular pulse width (ms)	0.98 (fixed)
Ventricular polarity pacing	BIPOLAR (fixed)

The device enters the "MRI Modality" with the programmed parameters when the follow-up visit is closed using the programmer; once the MRI examination is finished, the MRI Modality can be exited automatically when the set maximum time expires, or manually by the user.

For the operative details refer to Section 11.

### 9.1 Medico MRI pacing system integrity verification

Before enabling the MRI Modality, the Medico programmer automatically checks for possible problems with the pacemaker or the leads which may compromise patient safety.

**The only two conditions which stop the "MRI Modality" feature enabling are:**

- **Insufficient battery charge condition, in the event that the device is at RRT (Recommended Replacement Time).**
- **Status of the leads (neither UNKNOWN nor BIP->UNI).**

In these cases, the software reports with a ban symbol the condition of risk, and prevents the enabling of the "MRI Modality".

The software also checks the following conditions, and reports with an alert symbol whenever they are not fulfilled:

- High battery impedance - The software reports with a warning the risk state, if the measured battery impedance is greater than 5000  $\Omega$ .
- Lead impedance out of range - The software reports with a warning the risk state, if any of the bipolar atrial or ventricular lead impedance measurements is greater than 2000  $\Omega$ , is less than 200  $\Omega$ , or is unavailable.
- High ventricular pacing capture threshold – The software reports with a warning the risk state, if the ventricular capture threshold is  $\geq 2.0$  V at a pulse width of 0.43 ms.

14

- Leads implanted more than 6 weeks ago – The software reports with a warning the risk state, if the leads have been implanted less than 6 weeks ago.
- Bipolar pacing polarity for both chambers – The software reports with a warning the risk state, if the pacing lead configuration is unipolar.

**The physician has the responsibility to decide to ignore the warning messages and to proceed in any case with the “MRI Modality” enabling.**

## 9.2 “MRI Modality” feature activation

To activate the “MRI Modality” feature, proceed as follows:

1. Select “Programming” from the submenu on the right side of the home screen;
2. Select “MRI Modality” to open the “MRI Modality Checklist” (figure 1);
3. If the patient satisfies all the “MRI Modality Checklist” window listed conditions, answer “YES” at the question “Do you agree that the patient performs an MRI exam?”, on the lower side of the screen (figure 1);

**Note:** It is possible to print the “MRI Modality Checklist” window by hitting the “PRINT/SAVE” button, and choosing the “CURRENT PAGE” option;

4. Press the “Proceed” button, on the lower-right side of the screen (figure 1) to enter the “MRI Modality Configuration” window (figure 2).
5. Set the “**MRI Modality**” parameter, on the upper-left side of the “MRI Modality Configuration” window, up to “ON” (figure 2) and select the values of the underlying parameters according to the patient’s clinical needs, like the “Maximum Time at End of MRI Scan”, which defines the time-out of MRI Modality operation period. The physician has the responsibility to choose the most appropriate values.
6. Once the operative configuration is chosen, press the “Program” button in the lower-center side of the screen (figure 2). **This way, the “MRI Modality” is programmed and will start functioning immediately after the follow-up is finished.**
7. In order to assess the appropriateness of the MRI Modality, it is necessary to make an operative check of the programmed configuration pressing “START” button and then “STOP” button on the right side. This will ensure no diaphragmatic stimulation can be induced by the programmed energy.

Now the device is ready for an MRI scan.

Figure 1

Figure 2

## 10. DURING “MRI MODALITY” OPERATION

### 10.1 Device considerations

Suspension of counters and diagnostics – Device diagnostics and counters are suspended when “MRI Modality” is enabled.

Suspension of Magnetic Mode – The device does not start asynchronous, fixed-rate bradycardia pacing in the presence of a magnet when “MRI Modality” is enabled.

### 10.2 Radiology considerations during the MRI scan

Image artifact and distortion – Medico MRI leads have proved to grant a minimal image distortion for areas surrounding the implanted leads, when the device is out of the field of view. In the presence of the device near the field of view, a significant

image distortion could occur. When selecting the imaging parameters and field of view, it is important to consider image artifacts and distortion resulting from the presence of the device and the leads. It is necessary to consider these factors also when interpreting the MRI images.

**Patient monitoring** – A proper patient monitoring during the MRI procedure must be assured, including an accurate continuous monitoring of the patient's hemodynamic function.

The MRI environment could interfere with the patient monitoring system, so it is strongly recommended the use of more than one of the following systems:

- Electrocardiography
- Pulse oximetry
- Non-invasive blood pressure measurements

In addition, maintain visual and verbal contact with the patient.

**Indispensable preparation for patient rescue** – During MRI scan, keep an external defibrillator available.

**Note:** If, during the MRI scan, the patient's hemodynamic function is compromised, stop the MRI scan and take the proper measures to restore the patient's hemodynamic function.

## 11. AFTER THE MRI SCAN

After performing an MRI scan, it is important to fulfil the following operations:

- Disable the "MRI Modality" feature (see Section 11.1). Disabling the MRI Modality either manually or by time-out, does not restore the initial pacing energy. Instead, high energy pacing is performed as a safety measure until energy reprogramming is performed.
- Check the pacing capture threshold and reprogram the pacing energy to appropriate values.
- Check the battery status and the leads impedance
- Check the correct functioning of the whole system, with particular reference to pacing and sensing properties.
- If appropriate, calibrate again the TVI signal and activate the related diagnostics and features.

### 11.1 Restoring the device to the pre-MRI configuration

The device stores in memory the "MRI Modality" configuration programmed, until it is disabled either manual or by time-out.

Once the MRI scan has been completed, "MRI Modality" should be disabled using the programmer; if the time set with the parameter "Maximum time to MRI scan end" expires, the restoring of the operating function before the exam is automatic, **except for the pulse amplitude, which remains at the value set for the MRI scan.**

To disable the "MRI Modality" feature, perform the operations described in the Olympia programmer reference manual.

This restores the device parameters to the configuration preceding the "MRI Modality" activation, **except the ventricular output amplitude, that must be set to the most appropriate value after a threshold analysis is performed.**

**Note:** as mentioned above, at the end of the MRI scan it is strongly recommended to check again the pacing capture threshold, the correct sensing of the intrinsic signals, and the data related to the battery and the leads.

If the TVI feature is in use, a recalibration of the TVI signal is recommended before the corresponding features and diagnostics are re-activated.

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